

COOLINGBIS®

Monopolar electrosurgical electrode




01

COOLINGBIS® Features

COOLINGBIS® is a monopolar electro-surgical electrode intended for haemostatic sealing, coagulation and cut (models with blade) of soft tissues.

COOLINGBIS® uses radiofrequency (RF) energy and an **internally cooled electrode** to facilitate surgical sealing, increasing safety and drastically reducing intraoperative bleeding.



02

COOLINGBIS® Models and References

WORKING MODE

Coagulation and cut in
laparotomy (DUAL)

MODEL

Short electrode of 3 mm with blade
Short electrode of 5 mm with blade
Short electrode of 8 mm with blade

REFERENCE

BIS-3C11
BIS-5C11
BIS-8C11

Coagulation in laparotomy
(COAG)

Short electrode of 3 mm without blade
Short electrode of 5 mm without blade
Short electrode of 8 mm without blade

BIS-3C01
BIS-5C01
BIS-8C01

Coagulation and cut in
laparoscopy (DUAL)


Long electrode of 3 mm with blade
Long electrode of 5 mm with blade

BIS-3L11
BIS-5L11


Coagulation in laparoscopy
(COAG)

Long electrode of 3 mm without blade
Long electrode of 5 mm without blade


BIS-3L01
BIS-5L01



COOLINGBIS®
Ref. BIS-5L11



COOLINGBIS®
Ref. BIS-3C01



COOLINGBIS®
Ref. BIS-3L01

03

COOLINGBIS® Benefits



DUAL working mode



The DUAL working mode incorporates a blade to coagulate and cut the tissue, avoiding the need of using other dissecting devices.



Internal cooling system

The cooled saline flows inside the electrode without touching the tissue, which allows:

- > Preventing carbonized tissue from sticking to the electrode, producing a lower incidence of smoke.
- > Increasing the safety of the product. Some of the incidents caused by similar devices correspond to burns in the surgical area caused by heated irrigated saline.
- > Generating a lesion circumscribed under the electrode surface.





Lesion depths

The range of lesion depths that can be achieved with COOLINGBIS® depends on the **power level**, the **contact type** (tip only vs. lateral) and the **contact time**, **both superficial and deep lesions can be produced**.

Internal cooling allows enhancing lesions depths up to 6 mm and seal vessels up to 7 mm diameter (depending on models and power levels), without increasing the risk of thermal lesion to nearby structures. This allows for the creation of an additional ablation margin, which can **prevent local hepatic recurrence compared to conventional technologies, specially when the resection margin of healthy tissue that the surgeon left is limited** ¹.

Combinable

COOLINGBIS® allows the possibility of combining RF-assisted coagulation and cut with other devices (e.g. ultrasonic).

Faster

The design allows the rapid creation of coagulation lines by placing the whole electrode on the tissue. The electrode tip allows for easy sealing of blood vessels.

User-friendly

The curved shape of the electrode eases the movement and displacement along the working area, allowing easy access to areas that need to be coagulated and to the points of possible bleeding.

Coagulation efficiency

Without a doubt, COOLINGBIS® is currently **one of the most efficient coagulation devices on the market**. The high coagulation power drastically reduces blood loss and should therefore reduce the need for transfusion. Its high coagulation power is particularly useful in laparoscopic approaches where resources for effective and rapid haemostasis are more limited.

04

COOLINGBIS® Applications

COOLINGBIS® is intended to be used, but not limited to, in:

- > Liver
- > Kidney
- > Pancreas
- > Spleen²

Since it allows coagulation and cut of the tissue, it is especially intended for partial or total resection of these organs, i.e. for radiofrequency-assisted surgical resection.

COOLINGBIS® is especially useful in clinical cases that involve:

Livers with poor hemostasis:

- > Cirrhotic livers
- > Livers that underwent a lot of chemotherapy

Difficult to control bleeding:

- > Large vessels

Additional oncological ablation margins:

- > Very important in cases where it is necessary to prioritize minor resection

Pancreas:

- > To decrease pancreatic fistula in distal pancreas resections
-

05

COOLINGBIS® Clinical Outcomes In Liver

1 Reduction of complications:

Need of clamping manoeuvres is reduced³.

Intraoperative time is reduced³.


Plane of coagulative necrosis allows the optimum sealing of vessels and ducts and as a result may reduce the risk of biliary leakage^{4,5}.

2 Less need of transfusion:

Great power of coagulation which drastically reduces blood loss:

Results of the randomised clinical trial to evaluate the impact of COOLINGBIS® on intraoperative blood loss during liver resection (AGEMED 312/08 EC).⁶

The box diagram shows blood loss during transection (A) and blood loss per transection area (B) in the control group, in which haemostasis is obtained using conventional techniques, and in the MRFC (monopolar radiofrequency coagulation) group, in which haemostasis is obtained using COOLINGBIS®. The box represents the interquartile range containing 50% of the values. The whiskers are lines that extend from the box to the larger and smaller values. The line in the box indicates the median. Both results are statistically significant (*p<.05).




05

3 Reduction of local recurrence:


On the basis of the argument that remaining malignant cells in the hepatic remnant are responsible for tumour relapse, M. Villamonte, F. Burdío, E. Pueyo et. al. aimed to demonstrate that additional coagulation of the hepatic surface with an efficient RF-based device (COOLINGBIS®) not only successfully achieved hemostasis but also had a favorable effect on local recurrence⁷.

Figure shows the Kaplan-Meier curve of local hepatic recurrence-free survival in patients with liver tumors with distance from the tumour to resection margin <10 mm (Log-rank test $p=0.049$).

The estimated 1-, 3-, and 5-year local recurrence free survival of the control and COOLINGBIS® group patients were 93.5%, 86.0%, 81.0% and 98.8%, 97.2%, 91.9% respectively ($p=0.049$). The COOLINGBIS™ Group was significantly associated with reduced local recurrence⁷.




(A) The liver after removal of the sample; the coagulated tissue (<->) and the resection margin (*) can be seen. (B) The liver sample from the same patient; the resection margin and the thickness of the coagulated tissue (<->), which is in contact with the metastasis (->), can be seen. The * shows the correct position of the sample. (C) Histological section of the resection margins of the sample (use the asterisk for the correct position). The resection margin is marked with green ink. You can see the tumour (T) and the coagulated parenchyma (P) in contact with the margin. The coagulated tumour did not prevent the correct evaluation of the invasion of the margins. (D) CT scan of the same patient 56 months after liver resection, no signs of liver local recurrence are observed. The remaining ablated tissue at the margin can be noted (red arrow)¹.



05


④ Shorter postoperative recovery:

Significant reduction of **postoperative stay (days)**:



Results of consecutive patients in the COOLINGBIS® series in reference study 312/08/EC³. Patient n°3, after liver transection and in the same surgical act, underwent anterior resection of the rectum with hysterectomy and partial cystectomy due to invasive colon neoplasia of these organs. The patient was operated 11 days later as a result of a leakage of colorectal anastomosis from which he evolved favourably.

⑤ Transection plane is completely regenerated after procedure:



COOLINGBIS®

References

1. Quesada, R. et al. The impact of radiofrequency-assisted transection on local hepatic recurrence after resection of colorectal liver metastases. *Surg. Oncol.* 26, 229–235 (2017).
2. Quesada, R. et al. Laparoscopic partial splenectomy for giant cyst using a radiofrequency-assisted device: a case report. *Surg. Case Reports* 2, 2–5 (2016).
3. Burdío, F. et al. A new single-instrument technique for parenchyma division and hemostasis in liver resection: A clinical feasibility study. *Am. J. Surg.* 200, e75–e80 (2010).
4. Stavrou, G. A., Donati, M., Fruehauf, N. R., Flemming, P. & Oldhafer, K. J. Liver resection using heat coagulative necrosis: Indications and limits of a new method. *ANZ J. Surg.* 79, 624–628 (2009).
5. Weber, J. C. et al. New technique for liver resection using heat coagulative necrosis. *Ann. Surg.* 236, 560–563 (2002).
6. Quesada, R. et al. Impact of monopolar radiofrequency coagulation on intraoperative blood loss during liver resection: a prospective randomised controlled trial. *Int. J. Hyperth.* 33, 135–141 (2017).
7. Villamonte, M. et al. The impact of additional margin coagulation with radiofrequency in liver resections with subcentimetric margin: can we improve the oncological results? A propensity score matching study. *Eur. J. Surg. Oncol.* 25–29 (2021).

COOLINGBIS®



Vecmedical Spain S.L.
Pol. Ind. Can Milans, nave 9
Montcada i Reixac, Barcelona 08110

Tel. +34 93 572 50 15
vecmedical@vecmedical.com
www.vecmedical.com