

Textile Implantate Textile implants

SERATEX® MR

PP

POLYPROPYLEN
POLYPROPYLENE



gefärbt / monofil / nicht resorbierbar
dyed / monofilament / non-absorbable

SERATEX® PA / SERATEX® PA MR

PP/PGACL

POLYPROPYLEN /
POLYGLYKOLSÄURE-CAPROLACTON
POLYPROPYLENE /
POLYGLYCOLIC ACID-CAPROLACTONE



gefärbt / monofil / teilresorbierbar
Der PGACL-Anteil resorbiert innerhalb 90-120 Tage.

dyed / monofilament / partially-absorbable
The PGACL-part absorbs within 90-120 days.

MR = MRT sichtbar
= MRI visible

Symbol
Symbol

Bezeichnung
Designation


Verpackungs-Einheit
Units per box








Porengröße
Pore size








Abmessung
Size

Bänderanzahl
Number of tapes

Bestell-Nummer
Code number

PP		nicht resorbierbar / non-absorbable				
	SlimSling® SERATEX® L MR	3	2 x 1,5 mm	0,3 x 90 cm mit Nadel HRX-27	-	STL00401

PP / PGACL		teilresorbierbar / partially-absorbable				
	SERATEX® A PA MR	3	2 x 4 mm	2 x 19 cm	-	STA01401
	SERATEX® A2 PA MR	3	2 x 4 mm	50 x 1,5 cm	-	STA01402
	SERATEX® A4 PA MR	3	2 x 4 mm	55 x 1,5 cm	-	STA01404
	SERATEX® B1 PA	3	2 x 4 mm	2 x 15 cm	2 flex	STB01001
	SERATEX® B2 PA MR	3	2 x 4 mm	2 x 18 cm	-	STB01402
	SERATEX® B3 PA MR	3	2 x 4 mm	4,5 x 26 cm	-	STB01403
	SERATEX® B4 PA MR	3	2 x 4 mm	3,5 x 28 cm	-	STB01404


Symbol Symbol	Bezeichnung Designation	Verpackungs-Einheit Units per box	Porengröße Pore size	Abmessung Size	Bänderanzahl Number of tapes	Bestell-Nummer Code number
PP / PGACL		teilresorbierbar / partially-absorbable				
	SERATEX® E1 PA	3	2 x 3 mm	5 x 29 cm	-	STE01001
	SERATEX® E2 PA	3	2 x 3 mm	5 x 22,5 cm	-	STE01002
	SERATEX® E3 PA	3	2 x 3 mm	5 x 18 cm	-	STE01003
	SERATEX® E8 PA	3	2 x 3 mm	5 x 22,5 cm	-	STE01008
	SERATEX® E9 PA	3	2 x 3 mm	8 x 22,5 cm	-	STE01009
	SERATEX® E10 PA	3	2 x 3 mm	16 x 18,5 cm	-	STE01010
	SERATEX® F1 PA MR	3	2 x 4 mm	4,5 x 28 cm	4 flex	STF01401


SERATEX® PA

LIEFERPROGRAMM /
DELIVERY PROGRAMME





SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8-12
95119 Naila/Germany


 + 49 9282 937-0


 + 49 9282 937-9369

Export Department:

 + 49 9282 937-230

 + 49 9282 937-9785

 info@serag-wiessner.de

 www.serag-wiessner.de